

Introducing EXCIDIO®

A Nitrous Oxide Conversion Unit

Designed to support more sustainable clinical practice.






ENTONOX[®] in clinical care

ENTONOX remains an essential and widely used analgesic across:

- Maternity
- Emergency care
- Dental and minor medical procedures

Safe, effective, and trusted in everyday practice.



Supporting evolving priorities

Healthcare providers are increasingly focused on:

- Reducing environmental impact of anaesthetic gases
- Improving workplace air quality

EXCIDIO has been developed to support these goals without changing clinical practice.






How EXCIDIO works in practice

- Converts over 95% of exhaled nitrous oxide into oxygen and nitrogen
- Set-up in under 20 minutes
- Simple, single-button operation
- Mobile, quiet, standalone unit

Integrates into existing workflows with minimal disruption.



Practical impact

- Supports reduction in NHS Scope 1 nitrous oxide emissions
- Helps improve working environments for staff
- Enables more sustainable use of ENTONOX

Maintaining the care patients rely on, while addressing environmental priorities.

If you would like to learn more about EXCIDIO contact us at medicalbda@boc.com

